

CDKN2B/p15 INK4b Polyclonal Antibody FITC Conjugated

Catalog No: #C92160FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDKN2B/p15 INK4b Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Rat
Immunogen Type	65-138/138
Conjugates	FITC
Target Name	CDKN2B/p15 INK4b
Other Names	CDKN2B/p15 INK4b
Accession No.	Uniprot:42772Gene ID:1030
Uniprot	42772
GeneID	1030
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only