

FAT Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C92161PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FAT Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Rat
Immunogen Type	1401-1500/4588
Conjugates	PE-Cy5
Target Name	FAT
Other Names	FAT
Accession No.	Uniprot:Q14517Gene ID:2195
Uniprot	Q14517
GeneID	2195
Calculated MW	506
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only