

VE Cadherin Polyclonal Antibody FITC Conjugated

Catalog No: #C92164FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VE Cadherin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	251-320/784
Conjugates	FITC
Target Name	VE Cadherin
Other Names	VE Cadherin
Accession No.	Uniprot:P33151Gene ID:1003
Uniprot	P33151
GeneID	1003
Calculated MW	81
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only