

## EIF1AY Polyclonal Antibody FITC Conjugated

Catalog No: #C92165FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EIF1AY Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms RatCow
Immunogen Type	21-120/144
Conjugates	FITC
Target Name	EIF1AY
Other Names	EIF1AY
Accession No.	Uniprot:O14602Gene ID:9086
Uniprot	O14602
GeneID	9086
Calculated MW	16
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only