

SCRG1 Polyclonal Antibody FITC Conjugated

Catalog No: #C92166FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SCRG1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms RatPig
Immunogen Type	31-98/98
Conjugates	FITC
Target Name	SCRG1
Other Names	SCRG1
Accession No.	Uniprot:O75711Gene ID:11341
Uniprot	O75711
GeneID	11341
Calculated MW	11
Formulation	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only