

testis antigen 1B Polyclonal Antibody Biotin Conjugated

Catalog No: #C92169Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	testis antigen 1B Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu an
Immunogen Type	75-180/180
Conjugates	Bio
Target Name	testis antigen 1B
Other Names	testis antigen 1B
Accession No.	Uniprot:P78358Gene ID:1485
Uniprot	P78358
GeneID	1485
Calculated MW	18
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only