

# CSFV Serine protease NS3 Polyclonal Antibody HRP Conjugated



Catalog No: #C92176HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CSFV Serine protease NS3 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	CSFV
Immunogen Type	1951-2050/3898
Conjugates	HRP
Target Name	CSFV Serine protease NS3
Other Names	CSFV Serine protease NS3
Accession No.	Uniprot:P19712Gene ID:920146
Uniprot	P19712
GeneID	920146
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only