

EGF Polyclonal Antibody HRP Conjugated

Catalog No: #C92183HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGF Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC
Species Reactivity	Hu
Immunogen Type	1-53/53 (971-1023/1207)
Conjugates	HRP
Target Name	EGF
Other Names	EGF
Accession No.	Uniprot:P01133Gene ID:1950
Uniprot	P01133
GeneID	1950
Calculated MW	5.8/133
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-F=1:50-200

Note: This product is for in vitro research use only