

IL-6 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92186Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-6 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Pig Cow Hor Sheep
Immunogen Type	51-150/212
Conjugates	Bio
Target Name	IL-6
Other Names	IL-6
Accession No.	Uniprot:P05231 Gene ID:3569
Uniprot	P05231
GeneID	3569
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only