

## IGF I Polyclonal Antibody APC Conjugated

Catalog No: #C92187APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IGF I Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms RatDog
Immunogen Type	108-195/195
Conjugates	APC
Target Name	IGF I
Other Names	IGF I
Accession No.	Uniprot:P05019Gene ID:3479
Uniprot	P05019
GeneID	3479
Calculated MW	7.7/21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only