

TIMP-1 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C92193PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TIMP-1 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	131-207/207
Conjugates	PE-Cy5
Target Name	TIMP-1
Other Names	TIMP-1
Accession No.	Uniprot:P01033Gene ID:7076
Uniprot	P01033
GeneID	7076
Calculated MW	21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only