

ICAM1 Polyclonal Antibody HRP Conjugated

Catalog No: #C92198HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ICAM1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rab
Immunogen Type	451-505/532
Conjugates	HRP
Target Name	ICAM1
Other Names	ICAM1
Accession No.	Uniprot:P05362Gene ID:3383
Uniprot	P05362
GeneID	3383
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Note: This product is for in vitro research use only