

Phospho-Mst2 (Thr180) Polyclonal Antibody FITC Conjugated

Catalog No: #C92205FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-Mst2 (Thr180) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms RatRab
Immunogen Type	RN(p-T)VI
Conjugates	FITC
Target Name	Phospho-Mst2 (Thr180)
Other Names	Phospho-Mst2 (Thr180)
Accession No.	Uniprot:Q13188 Gene ID:6788
Uniprot	Q13188
GeneID	6788
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only