

NR2C2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92206Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NR2C2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-596/596
Conjugates	Bio
Target Name	NR2C2
Other Names	NR2C2
Accession No.	Uniprot:P49116Gene ID:7182
Uniprot	P49116
GeneID	7182
Calculated MW	65-67
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only