

## AMH Polyclonal Antibody Biotin Conjugated

Catalog No: #C92210Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AMH Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Ra
Immunogen Type	501-560/560
Conjugates	Bio
Target Name	AMH
Other Names	AMH
Accession No.	Uniprot:P03971Gene ID:268
Uniprot	P03971
GeneID	268
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only