

Macrophage mannose receptor 1 Polyclonal Antibody APC Conjugated



Catalog No: #C92218APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage mannose receptor 1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/1456
Conjugates	APC
Target Name	Macrophage mannose receptor 1
Other Names	Macrophage mannose receptor 1
Accession No.	Uniprot:P22897Gene ID:4360
Uniprot	P22897
GeneID	4360
Calculated MW	160
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only