

# Von Willebrand Factor Polyclonal Antibody Biotin Conjugated



Catalog No: #C92226Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Von Willebrand Factor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF
Species Reactivity	Hu Ms Ra
Immunogen Type	351-450/2813
Conjugates	Bio
Target Name	Von Willebrand Factor
Other Names	Von Willebrand Factor
Accession No.	Uniprot:P04275Gene ID:7450
Uniprot	P04275
GeneID	7450
Calculated MW	81
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only