

CD19 Polyclonal Antibody APC Conjugated

Catalog No: #C92227APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD19 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/556
Conjugates	APC
Target Name	CD19
Other Names	CD19
Accession No.	Uniprot:P15391Gene ID:930
Uniprot	P15391
GeneID	930
Calculated MW	59
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only