

CK15 Polyclonal Antibody PE Conjugated

Catalog No: #C92228PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CK15 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	351-456/456
Conjugates	PE
Target Name	CK15
Other Names	CK15
Accession No.	Uniprot:P19012Gene ID:3866
Uniprot	P19012
GeneID	3866
Calculated MW	49
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only