

CA153 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92229Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------------------------|
| Product Name | CA153 Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB |
| Species Reactivity | Hu |
| Immunogen Type | 941-1040/1255 |
| Conjugates | Bio |
| Target Name | CA153 |
| Other Names | CA153 |
| Accession No. | Uniprot:P15941Gene ID:4582 |
| Uniprot | P15941 |
| GeneID | 4582 |
| Calculated MW | 138 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200

Note: This product is for in vitro research use only