

# KLK3/Prostate Specific Antigen Polyclonal Antibody Cy3 Conjugated



Catalog No: #C92240Cy3

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KLK3/Prostate Specific Antigen Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-261/261
Conjugates	Cy3
Target Name	KLK3/Prostate Specific Antigen
Other Names	KLK3/Prostate Specific Antigen
Accession No.	Uniprot:P07288Gene ID:354
Uniprot	P07288
GeneID	354
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only