

LC3B Polyclonal Antibody FITC Conjugated

Catalog No: #C92253FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | LC3B Polyclonal Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1-50/125 |
| Conjugates | FITC |
| Target Name | LC3B |
| Other Names | LC3B |
| Accession No. | Uniprot:Q9GZQ8Gene ID:81631 |
| Uniprot | Q9GZQ8 |
| GeneID | 81631 |
| Calculated MW | 14 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only