

GOLPH2 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C92259Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GOLPH2 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/401
Conjugates	Cy5
Target Name	GOLPH2
Other Names	GOLPH2
Accession No.	Uniprot:Q8NBJ4Gene ID:51280
Uniprot	Q8NBJ4
GeneID	51280
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only