

CD24 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92262Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD24 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms
Immunogen Type	27-59/80
Conjugates	Bio
Target Name	CD24
Other Names	CD24
Accession No.	Uniprot:P25063Gene ID:100133941
Uniprot	P25063
GeneID	100133941
Calculated MW	44263
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only