

CD66a Polyclonal Antibody APC Conjugated

Catalog No: #C92265APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD66a Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	241-310/526
Conjugates	APC
Target Name	CD66a
Other Names	CD66a
Accession No.	Uniprot:P13688Gene ID:634
Uniprot	P13688
GeneID	634
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only