

KEAP1 Polyclonal Antibody PE Conjugated

Catalog No: #C92267PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KEAP1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/624
Conjugates	PE
Target Name	KEAP1
Other Names	KEAP1
Accession No.	Uniprot:Q14145Gene ID:9817
Uniprot	Q14145
GeneID	9817
Calculated MW	69
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only