

## MPO Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C92274PE-Cy5.5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MPO Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	678-745/745
Conjugates	PE-Cy5.5
Target Name	MPO
Other Names	MPO
Accession No.	Uniprot:P05164Gene ID:4353
Uniprot	P05164
GeneID	4353
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only