

CCK8 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92280Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CCK8 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200
Species Reactivity	Hu Ms Rat
Immunogen Type	96-103/115
Conjugates	Bio
Target Name	CCK8
Other Names	CCK8
Accession No.	Uniprot:P06307Gene ID:885
Uniprot	P06307
GeneID	885
Calculated MW	12.6
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P IHC-F ICC IF

Note: This product is for in vitro research use only