

Cystatin A Polyclonal Antibody PE Conjugated

Catalog No: #C92309PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Dog Cow
Immunogen Type	1-80/98
Conjugates	PE
Target Name	Cystatin A
Other Names	Cystatin A
Accession No.	Uniprot:P01040Gene ID:1475
Uniprot	P01040
GeneID	1475
Calculated MW	11
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only