

TFEB Polyclonal Antibody FITC Conjugated

Catalog No: #C92310FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TFEB Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/476
Conjugates	FITC
Target Name	TFEB
Other Names	TFEB
Accession No.	Uniprot:P19484Gene ID:7942
Uniprot	P19484
GeneID	7942
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only