

Alpha Dystroglycan Polyclonal Antibody Biotin Conjugated



Catalog No: #C92314Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha Dystroglycan Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	451-550/895
Conjugates	Bio
Target Name	Alpha Dystroglycan
Other Names	Alpha Dystroglycan
Accession No.	Uniprot:Q14118Gene ID:1605
Uniprot	Q14118
GeneID	1605
Calculated MW	98
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only