

# phospho-DOK1 (Tyr362) Polyclonal Antibody PE Conjugated



Catalog No: #C92326PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-DOK1 (Tyr362) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PI(p-Y)DE
Conjugates	PE
Target Name	phospho-DOK1 (Tyr362)
Other Names	phospho-DOK1 (Tyr362)
Accession No.	Uniprot:Q99704Gene ID:1796
Uniprot	Q99704
GeneID	1796
Calculated MW	62
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only