

Phospho-mTOR (Thr2446) Polyclonal Antibody PE Conjugated

Catalog No: #C92329PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-mTOR (Thr2446) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TR(p-T)DS
Conjugates	PE
Target Name	Phospho-mTOR (Thr2446)
Other Names	Phospho-mTOR (Thr2446)
Accession No.	Uniprot:P42345Gene ID:2475
Uniprot	P42345
GeneID	2475
Calculated MW	289
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only