

# phospho-SLK (Tyr530) Polyclonal Antibody PE Conjugated

Catalog No: #C92330PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-SLK (Tyr530) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PQ(p-Y)QP
Conjugates	PE
Target Name	phospho-SLK (Tyr530)
Other Names	phospho-SLK (Tyr530)
Accession No.	Uniprot:P06241Gene ID:2534
Uniprot	P06241
GeneID	2534
Calculated MW	59
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only