

# phospho-GFAP (Ser8) Polyclonal Antibody Biotin Conjugated

Catalog No: #C92331Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | phospho-GFAP (Ser8) Polyclonal Antibody Biotin Conjugated                       |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | WB IHC ICC IF   |
| Species Reactivity | Hu Ms RatPig  |
| Immunogen Type     | IT(p-S)A  |
| Conjugates         | Bio   |
| Target Name        | phospho-GFAP (Ser8)   |
| Other Names        | phospho-GFAP (Ser8)   |
| Accession No.      | Uniprot:P14136Gene ID:2670  |
| Uniprot            | P14136  |
| GeneID             | 2670  |
| Calculated MW      | 48  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only