

phospho-GSK-3 Beta (Ser21/Ser29) Polyclonal Antibody FITC Conjugated

Catalog No: #C92332FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-GSK-3 Beta (Ser21/Ser29) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	QP(p-S)AFGSMKV(p-S)RD
Conjugates	FITC
Target Name	phospho-GSK-3 Beta (Ser21/Ser29)
Other Names	phospho-GSK-3 Beta (Ser21/Ser29)
Accession No.	Uniprot:P49841Gene ID:2932
Uniprot	P49841
GeneID	2932
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only