

Phospho-c-Kit (Tyr721) Polyclonal Antibody FITC Conjugated

Catalog No: #C92338FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-c-Kit (Tyr721) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	NE(p-Y)MD
Conjugates	FITC
Target Name	Phospho-c-Kit (Tyr721)
Other Names	Phospho-c-Kit (Tyr721)
Accession No.	Uniprot:P10721Gene ID:3815
Uniprot	P10721
GeneID	3815
Calculated MW	105
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only