

PPAP2A Polyclonal Antibody PE Conjugated

Catalog No: #C92342PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAP2A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-130/284
Conjugates	PE
Target Name	PPAP2A
Other Names	PPAP2A
Accession No.	Uniprot:O14494Gene ID:8611
Uniprot	O14494
GeneID	8611
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only