

# phospho-IGF1R (Tyr1165 / Tyr1166) Polyclonal Antibody PE Conjugated

Catalog No: #C92343PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-IGF1R (Tyr1165 / Tyr1166) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TD(p-Y)(p-Y)RK
Conjugates	PE
Target Name	phospho-IGF1R (Tyr1165 / Tyr1166)
Other Names	phospho-IGF1R (Tyr1165 / Tyr1166)
Accession No.	Uniprot:P08069Gene ID:3480
Uniprot	P08069
GeneID	3480
Calculated MW	150
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only