

phospho-Ron (Tyr1238) Polyclonal Antibody PE Conjugated

Catalog No: #C92349PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-Ron (Tyr1238) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	E(p-Y)YS
Conjugates	PE
Target Name	phospho-Ron (Tyr1238)
Other Names	phospho-Ron (Tyr1238)
Accession No.	Uniprot:Q04912Gene ID:4486
Uniprot	Q04912
GeneID	4486
Calculated MW	30/119/150
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only