

phospho-PIK3R1 (Tyr368) Polyclonal Antibody FITC Conjugated



Catalog No: #C92352FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-PIK3R1 (Tyr368) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	GD(p-Y)TL
Conjugates	FITC
Target Name	phospho-PIK3R1 (Tyr368)
Other Names	phospho-PIK3R1 (Tyr368)
Accession No.	Uniprot:P27986Gene ID:5295
Uniprot	P27986
GeneID	5295
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only