

phospho-AMPK alpha-1 (Thr183) Polyclonal Antibody FITC Conjugated



Catalog No: #C92355FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-AMPK alpha-1 (Thr183) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	LR(p-T)SC
Conjugates	FITC
Target Name	phospho-AMPK alpha-1 (Thr183)
Other Names	phospho-AMPK alpha-1 (Thr183)
Accession No.	Uniprot:Q13131Gene ID:5562
Uniprot	Q13131
GeneID	5562
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only