

# phospho-PI3KCA (Tyr317) Polyclonal Antibody HRP Conjugated



Catalog No: #C92357HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-PI3KCA (Tyr317) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	RT(p-S)CG
Conjugates	HRP
Target Name	phospho-PI3KCA (Tyr317)
Other Names	phospho-PI3KCA (Tyr317)
Accession No.	Uniprot:P42336Gene ID:5290
Uniprot	P42336
GeneID	5290
Calculated MW	117
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only