

# phospho-AMPK alpha 2 (Ser173) Polyclonal Antibody FITC Conjugated



Catalog No: #C92358FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-AMPK alpha 2 (Ser173) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	RT(p-S)CG
Conjugates	FITC
Target Name	phospho-AMPK alpha 2 (Ser173)
Other Names	phospho-AMPK alpha 2 (Ser173)
Accession No.	Uniprot:P54646Gene ID:5563
Uniprot	P54646
GeneID	5563
Calculated MW	61
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only