

phospho-c-Raf (Ser43) Polyclonal Antibody FITC Conjugated

Catalog No: #C92364FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-c-Raf (Ser43) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	RA(p-S)DD
Conjugates	FITC
Target Name	phospho-c-Raf (Ser43)
Other Names	phospho-c-Raf (Ser43)
Accession No.	Uniprot:P04049Gene ID:5894
Uniprot	P04049
GeneID	5894
Calculated MW	71
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only