

phospho-RPS6KB1 (Ser427) Polyclonal Antibody APC Conjugated



Catalog No: #C92367APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-RPS6KB1 (Ser427) Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	KF(p-S)FE
Conjugates	APC
Target Name	phospho-RPS6KB1 (Ser427)
Other Names	phospho-RPS6KB1 (Ser427)
Accession No.	Uniprot:P23443Gene ID:6198
Uniprot	P23443
GeneID	6198
Calculated MW	58
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only