

ADAM12 Polyclonal Antibody FITC Conjugated

Catalog No: #C92393FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADAM12 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/909
Conjugates	FITC
Target Name	ADAM12
Other Names	ADAM12
Accession No.	Uniprot:O43184Gene ID:8038
Uniprot	O43184
GeneID	8038
Calculated MW	77/100
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=3ug/Test IF=1:50-200

Note: This product is for in vitro research use only