

ADAMTS13 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92397Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADAMTS13 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/1427
Conjugates	Bio
Target Name	ADAMTS13
Other Names	ADAMTS13
Accession No.	Uniprot:Q76LX8Gene ID:11093
Uniprot	Q76LX8
GeneID	11093
Calculated MW	145
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only