

Ferritin Heavy Chain/FTH1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92406Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ferritin Heavy Chain/FTH1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/183
Conjugates	Bio
Target Name	Ferritin Heavy Chain/FTH1
Other Names	Ferritin Heavy Chain/FTH1
Accession No.	Uniprot:P02794Gene ID:2495
Uniprot	P02794
GeneID	2495
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only