

AIM2 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C92415PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AIM2 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	255-354/354
Conjugates	PE-Cy5
Target Name	AIM2
Other Names	AIM2
Accession No.	Uniprot:Q91VJ1Gene ID:383619
Uniprot	Q91VJ1
GeneID	383619
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only